

PHARMACEUTICAL LABORATORY EQUIPMENT

Product Catalogue | Research | Quality Control | Teaching Laboratories



RESEARCH



QUALITY
CONTROL



TEACHING
LABORATORIES



PHARMACEUTICAL
SOLUTIONS

*Enabling
Science for
a Healthier
Tomorrow*



EQUIPMENT
FOR PHARMACY COLLEGES,
INSTITUTIONAL LABORATORIES,
RESEARCHERS AND



2026 EDITION | COMPLETE PRODUCT DIRECTORY & HISTORICAL PRICE LIST



centoptics.com

ambalasciencemarket.com



9812204116



centralworks17@gmail.com
kitchenindiav@gmail.com

PHARMACY LABORATORY EQUIPMENT DIRECTORY

For pharmacy colleges, R&D laboratories, quality-control teams and institutional science buyers.

Product codes, descriptions and historical prices are retained from the supplied 2020 CENTRAL catalogue.

IMPORTANT PRICE NOTICE

All prices shown in this edition are from the original 2020 list.

They are NOT confirmed current 2026 selling prices.

GST, delivery and availability must be verified before quotation.

PRODUCT ENQUIRIES

Phone: 9812204116

Email: centralworks17@gmail.com

Web: centoptics.com | ambalasciencemarket.com

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 001/A	ABBE'S REFRACTOMETER Measures R.I. 1.3 to 1.7; Accuracy 0.001 direct on scale & by estimation. Sugar percentage 0 to 95% & accuracy 1deg C by estimation. Measures dispersion & water circulation for constant temperature. Complete with test pieces, constant liquid thermometer with metallic jacket & Instruction manual.	Rs. 11,600.00
CW 001/B	HAND REFRACTOMETER (ERMA); 0-32%, 28-62%, 58-92% (Any one)	Rs. 3,200.00
CW 001/C	DIGITAL REFRACTOMETER "ERMA MAKE"	
PC001/C1	0-95% Brix With Refractive Index	Rs. 60,000.00
PC001/C2	0-280% Salinity & Refractive Index	P.O.R.
CW 002	ACTIVITY CAGE (Actophotometer) Consists of 6 built in photo sensor & 4 digit counter. Measures spontaneous & induced activity of rats with digital totalizer.	Rs. 10,800.00
PC 003	AMPOULE CLARITY TEST EQUIPMENT	Rs. 2,800.00
PC 004	AMPOULE FILLING AND SEALING DEVICE Syringe type motor driven giving ten ampoules per minutes. Capacity of filling adjustable from 1ml to 10ml by a cam lever arrangement. Automatic intake & outflow of fluid. All contact parts are of borosilicate glass. Consist of glass fusing burner for LPG with foot Below & Rubber tubing.	Rs. 8,800.00
PC 005	AMPOULE WASHING DEVICE Six jet nozzles & pan made of S.S. Spring loaded ampoule holding arrangement for various size of ampoules & tap to connect water from outside.	Rs. 3,600.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 006	ANALGESIOMETER	
PC 006/1	Radiant Heat type or tail Flick type	Rs. 6,400.00
PC 006/2	Eddy's Hot Plate	Rs. 10,400.00
PCN 01	ANDERSON PIPETTE (Borosilicate Glass)	Rs. 1,900.00
PC 007	ANTIBIOTIC ZONE READER (Manual)	Rs. 18,000.00
PCN 02	ARSENIC TEST APPARATUS	Rs. 300.00
PC 009	ASEPTIC CABINET (Wooden)	
PC 009/1	1.5' U.V. tube & one fluorescent tube	Rs. 9,600.00
PC 009/2	3' U.V. tube size & one fluorescent tube	Rs. 11,000.00
PC 010	AUTOCLAVE (Portable) 12"	
PC 010/1	Aluminium	Rs. 6,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 010/2	Stainless Steel	Rs. 9,000.00
PC 010/B	AUTOCLAVE (Vertical) Double Wall Inner Chamber 304 Grade S.S. & Outer M.S. Provided with Wing nut type locking Arrangement, Spring Loaded Safety Valve, Pressure Gauge, Pressure Release Valve, Water Level Indicator & drain Valve. Hydraulically Tested to withstand upto 40psi Pressure. Molded Jointless Silicon Rubber Gasket is Provided for leak Proofing.	
PC 010/3	250 x 450mm, 22 Ltr.	Rs. 40,000.00
PC 010/4	300 x 500mm, 40 Ltr.	Rs. 46,000.00
PC 010/5	350 x 550mm, 50 Ltr.	Rs. 56,000.00
PC 010/6	400 x 600mm, 75 Ltr.	Rs. 70,000.00
PC 010/20	450 x 600mm, 95 Ltr.	Rs. 80,000.00
PC 010/30	550 x 750mm, 187 Ltr. Optional Accessories:-	Rs. 120,000.00
PC 010/7	Automatic Low-Water Level Cut Off-Device.	Rs. 10,000.00
PC 010/8	Automatic Pressure Control Switch.	Rs. 10,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 010/9	Digital Temp. Indicator cum Controller	Rs. 8,600.00
PC 010/10	Radical lock	Rs. 9,000.00
PC 011	BALL MILL (Electrically Operated) Having steel balls of different sizes to be put in jar for mixing & grinding powder. (S.S. Jar)	
PC 011/1	Cap 1 Kg	Rs. 7,000.00
PC 011/2	Cap 2 Kg	Rs. 8,000.00
PCN 03	BOTTLE SEALING MACHINE (Hand Operated)	Rs. 19,600.00
PCN 04	BOTTLE WASHING MACHINE	
PCN 04/1	Table Top Model	Rs. 14,000.00
PCN 04/2	Floor Standing Model	Rs. 22,000.00
PCN 100	BROOKFIELD TYPE VISCOMETER	Rs. 260,000.00
PCN 05	BULK DENSITY APPARATUS	

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PCN 05/1	Motor 1/6HP, Cylinder cap. 50ml , Rpm Display : 3 digit LED	Rs. 11,400.00
PCN 05/2	With Digital Counter	Rs. 15,600.00
PCN 06	CANNULA (Glass)	
PCN 06/1	Symes Cannula	Rs. 40.00
PCN 06/2	Venus Cannula	Rs. 40.00
PCN 06/3	Lever Cannula	Rs. 40.00
PCN 150	CAPSULE COUNTER (Aluminium) 100 places	Rs. 2,400.00
PCN 07	CLAVENGER APPARATUS, Cap 1000ml (Borosilicate Glass)	Rs. 1,500.00
PC 013	POLE CLIMBING APPARATUS (COOK'S) Outer chamber made of thick plywood with observation window & inner made of acrylic sheet, with Brass C.P. bars at bottom. A shock for controlled duration of 200V AC 50Hz single phase 0.2 mA is built in & controlled by solid state timer.	Rs. 14,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 014	CAPSULE FILLING MACHINE (All Contact parts made of stainless steel.) A versatile lever operated semi automatic machine for various operations. (i) Plates for putting of capsules. (iii) Filling of powder by dispensing brush, pressing (ii) Removal of caps (iv) Ejection of 100 capsules of 0,00,1, 2, or 3 Any one size on order)	Rs. 28,000.00
PC 015	COLLAPSIBLE TUBE FILLING MACHINE (Hand Operated) Stroke type, barrel and piston arrangement with stroke setting for various capacity of filling tubes detachable nozzles provided for each separate capacity. All contact parts are S.S. Built on cast iron body. Cap. 5gm to 50gm	Rs. 9,200.00
PC 018	CONICAL PERCOLATOR (S. S. Jar on stand) 1000ml	Rs. 700.00
PC 019	ELECTRO CONVULSIOMETER Provides 50Hz stimulus current adjustable from 0.15 to 330mA in discrete steps for producing minimal and super maximal seizures required in study of anticonvulsant & epileptic drugs. Complete with accessories.	Rs. 9,200.00
PCN 08/1	DIGITAL POTENTIOMETER MAX Make	Rs. 8,200.00
PC 020	DISINTEGRATOR (Hopper Dia 2")	Rs. 16,400.00
PC 021	DISINTEGRATION TEST UNIT	
PC 021/1	Single Basket (Analog)	Rs. 7,200.00
PC 021/2	Double Basket (Analog)	Rs. 10,400.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 021/3	(Digital) Double Basket	Rs. 26,000.00
PC 021/4	(Digital) Single Basket	Rs. 22,000.00
PC 022	DISSOLUTION RATE TEST APPARATUS Consist of thermostatic control water bath which maintains temp. at 37DegC with + 0.5 Deg C variations. Provided with borosilicate vessel 1000 ml, variable speed motor, cylindrical wire mesh basket & shaft with blades of stainless steel.	
PC 022/1	Single Basket, (Analog)	Rs. 10,400.00
PC 022/2	Single Basket, (Digital rpm)	Rs. 14,600.00
PC 022/3	Single Basket, El Make (Model 912)	Rs. 44,000.00
PC 022/4	Six Test, With Digital timer, temp. & Rpm & Lift System	Rs. 120,000.00
PC 023	DISTILLATION PLANT (Manasty Type), Wall Hanging	
PC 023/1	Stainless Steel, Output approx. 3/4 Ltr./Hr.	Rs. 6,000.00
PC 023/2	Stainless Steel, Output approx. 5/6 Ltr./Hr.	Rs. 10,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 024	DISTILLATION APPARATUS (ALL GLASS) Automatic Electrically Heated All Glass Distillation Unit, NPL Design, Heater Embedded In Spiral Glass Tube Complete With Heavy Base, Rod, Clamps, Boards Etc.	
PC 024/1	Single Stage having flask of cap. 3 Litre	Rs. 12,000.00
PC 024/2	Single Stage having flask of cap. 5 Litre	Rs. 15,000.00
PC 024/3	Single Stage having flask of cap. 10 Litre	Rs. 22,000.00
PC 024/4	Double Stage having flask of cap. 3 Litre	Rs. 24,000.00
PC 024/5	Double Stage having flask of cap. 5 Litre	Rs. 30,000.00
PC 024/6	Double Stage having flask of cap. 10 Litre	Rs. 44,000.00
PC 024/10	Triple Stage having flask of cap. 3 Litre	Rs. 36,000.00
PC 024/11	Triple Stage having flask of cap. 5 Litre	Rs. 45,000.00
PC 024/12	Triple Stage having flask of cap. 10 Litre	Rs. 66,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 024/7	Automatic cut-off device for Single distillation	Rs. 3,600.00
PC 024/8	Automatic cut-off device for Double distillation	Rs. 4,400.00
PC 024/9	Automatic cut-off device for Triple distillation	Rs. 6,000.00
PC 024/20	Spare Flask with heater 3 Litre	Rs. 7,000.00
PC 024/30	Spare Flask with heater 5 Litre	Rs. 9,000.00
PC 024/40	Spare Flask with heater 10 Litre	Rs. 16,000.00
PC 240	DISTILLATION WATER STILL (BARNSTEAD TYPE) TABLE MODEL S.S	
PC 240/1	Output Capacity 4 to 5 Ltr/Hr. (Approx.)	Rs. 20,000.00
PC 240/2	Output Capacity 8 to 9 Ltr/Hr. (Approx.)	Rs. 40,000.00
PC 240/3	Output Capacity 15 to 16 Ltr/Hr. (Approx.)	Rs. 72,000.00
PCN 09	DOUBLE CONE BLENDER	

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PCN 09/1	5kg Cap.	Rs. 26,000.00
PCN 09/2	10kg Cap.	Rs. 42,000.00
PC 025	ELECTRONIC STIMULATOR (0- 30V)	Rs. 3,000.00
PCN 10	FILTER APPARATUS	
PCN 10/1	(Millipore Type) Consists of Membrane support, funnel, flask & clamp. Standard Joint fitting.	Rs. 3,600.00
PCN 10/2	(Simple Type) Consists of Porcelain Buchner Funnel fitted with Filter Flask by rubber bung.	Rs. 500.00
PCN 11	FLOURIMETER (Digital)	Rs. 72,000.00
PCN 12	FOLIN-WU TUBE	Rs. 70.00
PC 026	HAND GRINDING MILL	Rs. 2,000.00
PCN 13	GRANULATING SIEVE S.S. sieve fixed in a wooden frame	Rs. 550.00
PC 027	HISTAMINE CHAMBER (Complete with sphygmomanometer & Nebulizer)	Rs. 6,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 028	HOT AIR UNIVERSAL OVEN Double walled, Outer chamber made of M.S. Powder Coated. Temp. upto 250deg C. Thermostatic controlled, Heating Element on three sides. Thermometer is provided on top. The scale is calibrated in 5 deg C Steps Inner Chamber Size Capacity	
PC 028/1	Aluminium 12"x12"x12" 28Ltr	Rs. 7,000.00
PC 028/2	Stainless Steel 12"x12"x12" 28Ltr	Rs. 10,200.00
PC 028/3	Aluminium 14"x14"x14" 45Ltr	Rs. 8,000.00
PC 028/4	Stainless Steel 14"x14"x14" 45Ltr	Rs. 11,000.00
PC 028/5	Aluminium 18"x18"x18" 95Ltr	Rs. 12,400.00
PC 028/6	Stainless Steel 18"x18"x18" 95Ltr	Rs. 16,000.00
PC 028/11	Aluminium 24"x24"x24" 224Ltr	Rs. 20,000.00
PC 028/12	Stainless Steel 24"x24"x24" 224Ltr Optional:-	Rs. 28,000.00
PC 028/7	Air Circulating Fan	Rs. 3,600.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 028/8	Timer 0-24 hours	Rs. 3,800.00
PC 028/10	Digital Temp. Indicator-Cum-Controller	Rs. 8,600.00
PC 029	INCUBATOR BACTERIOLOGICAL Double walled, Outer chamber made of M.S. Powder Coated. Temp. from Ambient upto 80deg C. Thermostatic controlled, Heating Element on three sides. Thermometer is provided on top. The scale is calibrated in 5 deg C Inner Chamber Size Capacity	
PC 029/1	Aluminium 12"x12"x12" 28Ltr	Rs. 7,000.00
PC 029/2	Stainless Steel 12"x12"x12" 28Ltr	Rs. 10,200.00
PC 029/3	Aluminium 14"x14"x14" 45Ltr	Rs. 8,000.00
PC 029/4	Stainless Steel 14"x14"x14" 45Ltr	Rs. 11,000.00
PC 029/5	Aluminium 18"x18"x18" 95Ltr	Rs. 12,400.00
PC 029/6	Stainless Steel 18"x18"x18" 95Ltr	Rs. 16,000.00
PC 029/11	Aluminium 24"x24"x24" 224Ltr	Rs. 20,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 029/12	Stainless Steel 24"x24"x24" 224Ltr Optional Accessories Same as PC 028	Rs. 28,000.00
PC 030	LAMINAR AIR FLOW UNIT Structure made of thick ply board with sunmica from outside & S.S. Workbench. Clear acrylic side panels & front door.	
PC 030/A	HORIZONTAL UNIT	
PC 030/1	60x60x60cm	Rs. 58,000.00
PC 030/2	90x60x60cm	Rs. 70,000.00
PC 030/3	120x60x60cm	Rs. 80,000.00
PC 030/B	VERTICAL UNIT	
PC 030/4	60x60x60cm	Rs. 70,000.00
PC 030/5	90x60x60cm	Rs. 82,000.00
PC 030/6	120x60x60cm	Rs. 99,000.00
PCN 14	LEAK TEST APPARATUS 34000.00	

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 031	LIPSTICK MOULD	Rs. 1,300.00
PCN 15	LIQUID / BOTTLE FILLING MACHINE	Rs. 19,000.00
PCN 16	LUCAS MOIST CHAMBER	Rs. 1,300.00
PCN 17	MECHANICAL AGITATOR (Propeller Type)	Rs. 7,000.00
PC 032	MILLIPORE Type FILTER 3 grades (Dia 47 mm) different positions.	Rs. 3,000.00
PC 033	OINTMENT SLAB (made of acrylic sheet size 8" x8")	Rs. 200.00
PC 034	OINTMENT SPATULA OR DRUG SPATULA (Stainless Steel)	
PC 034/1	6" wide blade, stout, flexible, with wooden	Rs. 60.00
PC 034/2	8" wide blade, stout, flexible, with wooden	Rs. 70.00
PC 035	OPERATION TABLE (Small) 30x30cms. S.S. Top	Rs. 5,000.00
PC 036	PARIEMETER	Rs. 15,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PCN 18	PERMEABILITY CUP	Rs. 1,600.00
PCN 19	PHASE CONTRAST WITH BINOCULAR MICROSCOPE	Rs. 88,000.00
PCN 20	PLETHYSMOGRAPH	Rs. 1,500.00
PC 037	POLISHING PAN Inner Dia 12 x8" high with opening of 7" Dia neck. For Use With Coating Pan. Made Of M.s. & Provided With Inner Lining Of Superior Canvas Cloth.	Rs. 3,400.00
PC 039	RABBIT CAGE size 18" x 18" x14" With heavy gauge mesh mounted in a frame of M.S angle with door in front. 2" Deep tray of G.I sheet is provided which slides in an angle to the sides of the cage and is removable for cleaning.	Rs. 3,600.00
PC 040	RABBIT HOLDER (Size: 17" x 6" x 6") Box type Wooden glass front, top lid opening air window, pushing system & neck holding clamp provided.	Rs. 2,200.00
PCN 21	REFLUX FLASK WITH CONDENSER (B-24 joint)	
PCN 21/1	Cap. 250ml	Rs. 450.00
PCN 21/2	Cap. 500ml	Rs. 500.00
PCN 21/3	Cap. 1000ml	Rs. 700.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PCN 200	RHEOMETER	Rs. 300,000.00
PC 041	ROTA ROD	
PC 041/1	Four Compartment	Rs. 16,000.00
PC 041/2	Three Compartment	Rs. 14,000.00
PC 041/3	Two Compartment	Rs. 13,000.00
PC 016	SOXHLET Glass Parts Only (Borosilicate glass)	
PC 016/1	Flask Cap. 250ml & Soxhlet Capacity 100ml	Rs. 800.00
PC 016/2	Flask Cap. 500ml & Soxhlet Capacity 200ml	Rs. 1,200.00
PC 016/3	Flask Cap. 1000ml & Soxhlet Capacity 500ml	Rs. 2,200.00
PC 017	SOXHELT EXTRACTIONS UNITS (Without Glass Parts) Cap upto 500ml, One Regulator	
PC 017/1	3 Mantles (Mantle Type)	Rs. 4,800.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 017/2	6 Mantles (Mantle Type)	Rs. 6,400.00
PC 017/8	Extra Regulator (Each)	Rs. 700.00
PCN 22	SIEVE SHAKER	
PCN 22/1	Hand Operated (For 8" and 12" Dia Sieve)	Rs. 8,000.00
PCN 22/2	Electrical Table Model (Only For 8" Dia Sieve)	Rs. 26,000.00
PCN 22/3	Gyratory Motion (Heavy Pattern) For 8" & 12" Dia Sieve	Rs. 39,000.00
PCN 22/4	Optional Timer For Above	Rs. 2,200.00
PC 043	SMOKING STAND	Rs. 2,800.00
PC 044	SMOKING BURNER	Rs. 1,400.00
PCN 23	SONICATOR (Ultra Sonic Bath), Complete S.S & Digital Temp. Controlled	
PCN 23/2	Capacity 2.5Ltr Approx	Rs. 20,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PCN 23/3	Capacity 4Ltr Approx	Rs. 37,000.00
PCN 23/4	Capacity 6 Ltr Approx	Rs. 56,000.00
PCN 23/B	PROBE SONICATOR	Rs. 370,000.00
PCN 24	SPRAYER (Glass) With Bellow	Rs. 400.00
PCN 25	SPIROMETER	Rs. 12,000.00
PCN 26	STERILITY TEST APPARATUS	Rs. 44,000.00
PC 045	ORGAN BATH	
PC 045/1	Double Unit	Rs. 6,600.00
PC 045/2	Single Unit	Rs. 2,700.00
PC 046	SUPPOSITORY MOULD (Brass)	
PC 046/1	4 Hole, 1gm	Rs. 450.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 046/2	6 Hole, 1gm	Rs. 500.00
PC 046/3	4 Hole, 1.5gm	Rs. 650.00
PC 046/4	6 Hole, 1.5gm	Rs. 850.00
PC 046/5	12 Hole, 1gm	Rs. 1,050.00
PC 047	TABLET MAKING MACHINE (Punching Machine) SINGLE PUNCH for making Tablets up to 1000 mg as per size of die & punch 6, 9 & 12 mm (Any one).Automatic compression & ejection device with each revolution of the wheel drive. Built on heavy cast iron body with eccentric pressure system. Stainless steel hopper provided & connected by a Rubber tube. Easy setting for thickness & hardness of tablets with lock nuts. Die &	
PC 047/1	Hand Operated	Rs. 26,000.00
PC 047/2	Electrically Operated	Rs. 44,000.00
PC 048	TABLET COATING PAN 12" Dia stainless steel, the machine worked by oil immersed gearbox, with F.H Motor drive. Air blower is operated independently & air is thermostatically controlled to give temperature 30 to 110 degrees. Switches for motor, blower, heater and Thermostat are fixed on a panel with neon indicators. Also mains lead with 15 amps plug provided.	Rs. 28,000.00
PC 049	TABLET COUNTER Suitable for 100 tablet	Rs. 1,100.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 050	TABLET FRIABILITY TEST APPARATUS The drum rotates at a speed of 30rpm by s geared motor & is controlled by a timer.	
PC 050/1	Single Drum (Analog)	Rs. 8,000.00
PC 050/2	Double Drum (Analog)	Rs. 11,000.00
PC 050/3	Single Drum DIGITAL	Rs. 19,000.00
PC 050/4	Double Drum DIGITAL	Rs. 22,000.00
PC 051	TABLET HARDNESS TESTER	
PC 051/1	Monsanto type	Rs. 1,500.00
PC 051/2	Pfizer Type	Rs. 4,400.00
PC 052	TELE THERMOMETER	
PC 052/1	with Single Probe	Rs. 7,600.00
PC 052/2	with 6 probe	Rs. 9,000.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 053	TINCHER PRESS (1 Ltr.) Contact Parts made of S.S.	Rs. 2,600.00
PC 055	COLLAPSIBLE TUBE CRIMPING AND SEALING MACHINE	Rs. 2,500.00
PC 056	VOLUME RECORDER	Rs. 1,400.00
PC 058	WATER DEIONIZATION UNIT (DEIONISER) Consists of two plastic column charged with high exchange capacity Cation & Anion resins. Conductivity Meter Provided .	
PC 058/1	Capacity 50 Ltr/Hr	Rs. 20,000.00
PC 058/2	Capacity 100 Ltr/Hr	Rs. 40,000.00
PC 058/3	Capacity 200 Ltr/Hr	Rs. 28,000.00
PC 059	X- BLOCK (Open)	Rs. 200.00
PC 060	ELECTRICAL E-8 STUDENT KYMOGRAPH (Recording Drum)	Rs. 12,000.00
PC ACC	GENERAL ACCESSORIES FOR KYMOGRAPH	

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 060/1	Adjustable Stand	Rs. 440.00
PC 060/2	Drum Paper Clip	Rs. 150.00
PC 060/3	Myograph Lever	Rs. 400.00
PC 060/4	Frog Board	Rs. 400.00
PC 060/6	Frog Muscle Chamber	Rs. 1,100.00
PC 060/7	Frog Heart Chamber	Rs. 1,100.00
PC 060/8	Induction coil Dubois Reymond	Rs. 3,400.00
PC 060/9	Dubois key	Rs. 800.00
PC 060/10	Muscle Holder or Feemur Clamp	Rs. 400.00
PC 060/11	Simple Electrode	Rs. 150.00
PC 060/12	Tuning Fork N-100	Rs. 400.00

INSTRUMENTS & EQUIPMENT

Catalogue model / product specification / historical price (INR)

MODEL CODE	PRODUCT & SPECIFICATIONS	PRICE (INR)
PC 060/13	Glazed Paper Sheet (100 sheets)	Rs. 400.00
PC 060/14	Starling Heart Lever	Rs. 200.00
PC 060/15	Brodie Universal Lever	Rs. 400.00
PC 060/16	Marry Tambour, 28mm	Rs. 500.00
PC 060/18	Artery Forceps	Rs. 200.00
PC 060/19	Frogs Pithing Needle	Rs. 70.00
PCN 27	Muscle Electrode	Rs. 400.00
PC 060/24	Heart Clip, S.S.	Rs. 200.00
PC 060/27	Aerator	Rs. 450.00
PC 060/30	Aeration Tube Borosilicate Glass	Rs. 70.00